





NEW SPECULAR MICROSCOPE

EM-4000




Why **Endothelium** Should Be Measured (1)

EM-4000

Endothelial Cell Density change by Aging



Age	CD (cells/mm ²)
New born	3,000 - 4,000
10 - 19	2,900 - 3,500
20 - 29	2,600 - 3,400
30 - 39	2,400 - 3,200
40 - 49	2,300 - 3,100
50 - 59	2,100 - 2,900
60 - 69	2,000 - 2,800
70 - 79	1,800 - 2,600
80 - 89	1,500 - 2,300



Why Endothelium Should Be Measured (2)

EM-4000

Various Clinical application



Pre-surgery:

- Keratorefractive surgery
- Cataract Surgery
- Phakic IOL surgery
- Contact Lens Fitting

Cornea disease:

- Physical forces
- Cornea guttata(Fuchs Dystrophy)
- Polymorphism

Post surgery:

- Keratoplasty
- IOL, Catarract surgery

New EM-4000 Features

EM-4000

1. **Newly** well-designed shape
2. TOMEY standardized screen menu
w/ **larger (10.4")** TFT touch screen
3. Para-central fixation lamp
(**Wider capturing**)
4. **Various analysis** methods
(TOMEY Center Method)
5. **Built-in printer** with auto-cutter
6. Internal **database** (SD card)
7. **Easy** operation and **Fast** analysis



TOMEY SPECULAR Microscope
EM-4000


TOMEY Standardized screen with Full-auto operation

EM-4000



OA-2000

- **Operation with new GUI**
 - Same layout as OA-2000
 - 10.4 inch touch screen
 - Same designed buttons
- **Automatic measurement**
 - Touch screen controlled
 - Auto-alignment
 - Automatic light adjustment
 - Auto-shot



EM-4000 operation screen


Analysis Speed (Quick measurement)

EM-4000



1~2 SEC PER EYE!!

EM-4000 measurement mode



Analysis Speed

Company	EM-4000	EM-3000	T company	N company	C company
Auto Analysis	1-2 sec	8-10 sec	Checking	3-5 sec	Checking

AUTO Trace/Analysis Method EM-4000



QUICK & EASY
YET FULL ANALYZING CHOICE

STEP1



EM-4000 measurement mode

STEP2




Select measured data for analysis

Analyzing modes of measured data



Measured data (Auto) trace




Cell size




Cell shape

STEP3



Right/Left Eye mode

OR



Single Eye mode

NEW Analysis Method (Core method)


EM-4000

TOMEY Core Method

NEW additional analysis method, TOMEY Core method, is available by pointing a center of cell for the count.(similar to KONAN Center Method)



Analysis screen after Core Method





Core Method screen

Company	EM-4000	EM-3000	T company	N company	C company
Center Method	Available	NOT available	NOT available	Available	NOT available

TOMEY L-cont Method (User' s specified area analyzation)

EM-4000

The cell in the L-shaped-range is counted for analyzation!




Company	EM-4000	EM-3000	T company	N company	C company
Center Method	YES	YES	NO	NO	NO


DARK area analysis mode

EM-4000

FURTHER ANALYSIS




EM-4000 analysis screen



Dark area

CCT	494	508
CCT (US)	507	521
DA NUM	0	7
DAD	0	623
DA AVG	0	1606
DA SD	0	1596
DA CV	0	99
DA MAX	0	5153
DA MIN	0	299
DA Ratio 0.0	10.7	
(D) (A)		

Print: 2014/Jun/19 17:54
TOMEY CORP.



DARK Area analysis screen
with Dark area number, Density,
Dark area ratio and more at glance


With EM-4000 setting,
Dark area analysis can be added to
printed information.

Para-central & Near-central w/ Fixation light position

EM-4000

Measuring area **0.25 x 0.54mm** is achieved with central points as well as **8 Near-central points** & **6 Para-central points**

ing peripheral.



Company	EM-4000	EM-3000	T Company	N Company	C Company
Para-Central	Available	NOT available	Available	Available	NOT available

Internal Data Base

EM-4000

With Internal data base on SD card,
About **16,000 patients** in 32GByte for analysis can be stored.

Comparison check with pre & post surgery can be checked on EM-4000 unit.



Company	EM-\$000	EM-3000	T Company	N Company	C company
Internal DB	YES	NO	NO	NO	Checking

Built-in Printer for diagnostic information

EM-4000

- **Easy printing**


For measured data and image data.



YMD:2014/Jun/19 17:54:21
EM-4000 Ver.00

ID :12345678
Name:TOMEY TARO
Sex :Male

RIGHT LEFT
FIX C I C




	* [Auto]	[Auto]
NUM	310	204
CD	2759	2181
AVG	362	459
SD	149	402
CV	41	88
MAX	1312	5153
MIN	73	104
6A	42	39

New Functions compared with EM-3000

EM-4000

- **Built-in printer**
 - Printing of analysis values and **cell-images**.
- **Database**
 - **SD card as storage**
 - Storage of approximately **16,000 patients**.
- **Additional Analysis: Core Method**
 - **Manual analysis of cell density**
 - **KONAN center method**
- **Fixation lights for Para-central exams**
 - **Additional 8 fixation lights**.



Customization of data output

EM-4000



(Fig. 1)

Functions that may be assigned to a button are as follows.



Export



Print



Save



Save & Export



Print & Export



Print & Save



ALL (Export, Print, and Save)

Comparison With EM-3000 (Improvements)

EM-4000

	EM-3000	EM-4000
Measurement areas (Possibility)	<ul style="list-style-type: none"> Center Peripheral 6 point 	<ul style="list-style-type: none"> Center Peripheral 6 points Near central 8 points
Analysis method	<ul style="list-style-type: none"> Auto-trace L-count Dark area 	<ul style="list-style-type: none"> Auto-trace L-count (manual/auto) Dark area Core method
internal Database function	Not Supported	SD card (32 GB) (16,000 patients)
Patient Follow-up	only with TB-1000	YES R&R L&L view
Built in printer	No	Yes
Monitor Size	8.4 inch	10.4 inch
Monitor tilt	Not Supported	Yes
Connectors	LAN	LAN / USB /USB PictBridge

Comparison With Competitors

EM-4000

	TOMEY EM-4000	N company	T company	K company	C company
Measurement image by 1 shot	16	16	1	1	9
Measurement range	0.54×0.25	0.55×0.25	0.55×0.25	0.4×0.24	0.54 x 0.27
Measurement Point	Center1 + periphery 6 points + Para central 8points	Center1 + periphery 6 points + Para central 8points	Center3 + periphery 6 points + Para central 8points	Center1 + periphery 4 points	Center1 + periphery 6 points
Possibility for manual analysis	L-count / Core method/Trace method	1 kind of Manual Analysis	Not Supported	Center method/Flex center method/ Trace method	1 kind of Manual Analysis
Database function	Supported	Not Supported	Not Supported	Not Supported	Supported
Language	ENG/JAP	ENG/JAP	ENG/JAP	ENG/JAP	Eng./German/ Italian/Spanish/ Portuguese/Dutch
Panoramic Measurement	Not Supported	Not Supported	Supported	Not Supported	Not Supported
Display Size	10.4 inch	8.4 inch	10.4 inch	15 inch	10.4inch

EM-4000



*Thank you for
your attention!*